

staxi® Medical Chair - 600lb Capacity

Anti-microbial Cushions

For optimum hygiene and ease of maintenance.

Lifting Armrests

Facilitates easy and safe side entry and exit of occupants. Anti-Microbial.

Non Folding Frame

Strong, durable and safe.

Footrest

For occupant comfort and safety.

Turns within its own circumference



Fail-safe Brake System

For occupant safety and security.

Courtesy Rack

Ideal for storage of small bags, papers, magazines. Hanging hooks 25lbs weight capacity – ideal for handbags, hand luggage.

Large Luggage Area

For occupant convenience and safety.

High Quality Castors & Anti-tip Design

Easy to manoeuvre for a smooth and reliable ride.

STAXI Medical Chair – Technical Specification

Overall Dimensions				
Width*	Length	Seat to Floor	Back Height	Seat Width
27"(686mm)	41" (1042mm)	21" (534mm)	40.5" (1029mm)	20.5" (521mm)

Nesting Dimensions		
Nest	Nest with cushion	First chair
11" (280mm)	12" (305mm)	41" (1042mm)

Weight		
Chair weight	Shipping Weight	Occupancy Capacity
57lbs	80lbs	600lbs

Colors
Automotive Grade. Specially Formulated. Mar-resistant. Undercoat & Overcoat. Weather Resistant
Blue, Red, Silver and Diamond Black.

Chair Specifications	
Arm Rests	Anti Microbial. Steel Oval Tubing. Non-porous sleeve. Lift for Easy Transfer
Footrests	Steel. Non-Removable. Lifts for comfort
Brake	Automatic; Fail Safe Ratchet
Wheels	2 1/8 inch composite wheel. 8 inch diameter. Non-marring
Castors	5-inch. Urethane. Sealed bearing. Non-marring
Frame	Heavy duty; Formed oval tubing
Cushion	Antimicrobial treated upholstery.

STAXI Corporation recommend that each STAXI transportation chair is fitted with a STAXI seatbelt. For further details contact customer services on 1.866.969.1463 or email customerservice@staxi.com

The **staxi®** Medical comes with a **3-year parts warranty**.

* **Overall Overall Chair Width Including O2 Holder: 30" (762mm)**



www.staxi.com
customerservice@staxi.com
 1.866.969.1463

staxi®
 Progress Through Innovation™