



BODY FLUID MANAGEMENT SYSTEMS



KERMA MEDICAL IS AN AUTHORIZED DISTRIBUTOR OF HYGIE PRODUCTS IN THE USA



BEDPAN

BEDPAN HYGIENIC COVER® SA-HYGI-CLBI-U00

The Bedpan/Commode Hygienic cover® is the perfect solution to effectively manage and contain body wastes at point of care. Its unique super-absorbent pad included in each cover (liner) is designed to turn more than 700 ml of body fluids (urine, blood and feces) into gel within seconds.

The Hygienic Covers® are made from recycled material and the super-absorbent pad is made from 100% recycled material.



BEDPAN SUPPORT BA-BASS-HY21-000

Used with Hygienic covers®, the Hygie Bedpan Support® is the perfect solution to contain pathogens. It combines the characteristics of the orthopedic bedpan and the standard bedpan. Your staff will never have to handle stool again, clean and disinfect soiled equipment. Recommendation: one support per patient per stay.

- Increased stability
- Easy handling
- Easy storage (stackable)
- No sterilization required
- 100% recyclable



BEDPAN

BEDPAN PROKIT KT-PRO-HY21-B00

Combines in a value format, the bedpan support, Bedpan/Commode Hygienic Covers®. A complete system to manage body fluids in your institution, whether you need to isolate or standardize.

- 3-in-1 Kit
- No cleaning
- No spills
- No sterilization
- Easy to implement



PADS

SUPER-ABSORBENT PAD SA-TAMP-ABSO-020

The super-absorbent pad turns 600 ml of body fluids into gel within 30 seconds. It is particularly recommended for nosocomial infections (HAI's), reducing the stress associated with the cleaning and emptying of traditional bedpans and the risks of spills and splashes. Convenient and effective, it makes the collection, transportation and disposal of body fluids safer. Made from 100% recycled materials, it is used inside the Hygienic covers® for bedpan/commode, vomit and urinal supports.

- Absorbs 600 ml or more in 30 sec
- Simple and easy to use
- No Cleaning
- No spills



URINAL

URINAL HYGIENIC COVER® SA-HYGI-URIO-000

The Urinal Hygienic Covers® made from recycled material are the fast and efficient solution for disposing of urine. The superabsorbent pad included with the covers will eliminate inconveniences such as cleaning and splashing while reducing the risk of cross-contamination.

- Includes a super-absorbent pad
- No cross-contamination
- No Transportation
- Made from recycled materials
- Recommended for all units or all types of institutions (Hospitals or Long Term Care)



URINAL SUPPORT HY-SUPP-URI3-000

The Urinal Support can effectively manage and contain urine at the point of care. Designed for use with Urinal Hygienic Covers®.

- Easy grip
- Unisex
- Splatter-free
- Increased safety
- 100% Recyclable



URINAL

URINAL PROKIT KT-PRO-HY21-ULO

The Urinal ProKit combines in a single value format the Urinal Support and the Urinal Hygienic covers®. The complete system to manage urine effectively in your facility.

- No cross-contamination
- No Transportation
- Made from recycled materials
- Economic
- 3-in-1 pack



VOMIT

VOMIT HYGIENIC COVER® SA-HYGI-VOM0-000

Hygienic Covers® for vomit supports are made from recycled material. They help contain contaminants at the source and include a super-absorbent pad that solidifies vomit in seconds eliminating spills and splashes.

- Includes a super-absorbent pad
- No cross-contamination
- No Transportation
- Made from recycled materials
- Recommended for all units or all types of situations (emergencies, oncology, medicine, ambulances, etc.)



VOMIT SUPPORT HY-SUPP-VOM0-000

Designed to be used with Hygienic Covers® (vomit bag) and to capture vomit in all situations. The super-absorbent pad inserted in each Hygienic Cover® solidifies body fluids in seconds eliminating spills and splashes.

It can be hooked to the bed side bar or commode and the multiple handles around the support allow for a better grip.

- Contaminant Isolation
- Easy-grip
- Multi-functional (Vide stomie)
- Easy storage (stackable)
- 100% recyclable



VOMIT

GO VOM VOMIT BAG – 35 BAGS SA-HYGI-GOVO-000

The GO VOM is a portable vomit bag with a cardboard cone at the opening, an anti-reflux valve and an absorbent pad. The super-absorbent pad inside solidifies over 300 ml of body fluids within seconds. This bag is graduated and transparent to facilitate measurements.

- Ideal for emergencies (ER, ambulance)
- Pads are made of 100% recycled material.
- Fast and easy to implement
- Useful for home care and Long Term Care
- No cleaning



VOMIT SUPPORT KT-PRO-HY21-V01

This ProKit effectively manages vomiting in any circumstance. It combines in a single economical format the *Vomit Support® and Vomit Hygienic covers® (vomit bags)*.

- No cleaning
- No spills
- Ideal for emptying colostomy bags
- Convenient for chronic pulmonary obstruction
- Fast and easy to implement

